

ASPEN VISTA® MULTIPOST COLLAR

MULTIPOST CERVICAL COLLAR
MULTIPOST ZERVIKALORTHESE
CUELLO CERVICAL MULTIPOST
COLLARE CERVICALE MULTI-POSIZIONE
COLLIER CERVICAL À PLUSIEURS POSTES
MULTIPOST NECKBRACE
VÍCE-POLOHOVÝ KRČNÍ LÍMEC
MULTIPOST HALSHVIRVELKRAVE
ПОЛУКОМБИКО АУХЕНИКО КОЛАРО
MULTIPOST CERVICAL KRAAG
KOLNIERZ ORTOPEDYCZNY WIELOPODPOROWY
MULTIPOST KAULATUKI
MULTIPOST NACKKRAGE
МНОГОПОЗИЦИОННЫЙ ШЕЙНЫЙ КОРСЕТ
COLAR CERVICAL MULTIPOST
多重支撑颈托

ASPEN VISTA CERVICAL COLLAR FRONT

ASPEN VISTA MULTIPOST BACK PANEL

MEDICAL DEVICE

SINGLE PATIENT - MULTIPLE USE

CAUTION


EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Pat. aspemp.com/patents

Aspen Medical Products, LLC
6481 Oak Canyon, Irvine, CA 92618-5202 USA
aspemp.com • +1-800-295-2776 • +1-949-681-0200

©2021 Aspen Medical Products, LLC. All Rights Reserved.
AC0113D revised 04/2021

1

ENGLISH

1. PREPARATION

Pre-form the back panel so plastic side pieces are slightly curved.
NOTE: Ensure head is in neutral alignment throughout the entire application of this product.

2. PLACEMENT

The front panel should be at the lowest setting to begin. While pulling the sides of the front panel apart, position the bottom of the tracheal aperture at the sternal notch. Position the sides up and over the trapezius muscles. It is OK if the chin piece is not touching the chin.

3. SIZING

While holding the collar against the chest, pull out the dial to unlock. Turn the dial clockwise to raise chin piece. When the chin piece is supporting the chin, release the dial and it will self-lock.

4. MULTIPOST BACK PANEL

NOTE: Do not adjust on the patient. If height adjustment is necessary - remove from patient, adjust and reapply.

A. Push down on lock release and slide yellow button left to unlock.
B. Adjust height up or down if needed for best fit. To lock, slide yellow button right.
C. With the head in neutral alignment, press the back panel down onto the bed and side it behind the patient's neck. Ensure back panel has good contact with trapezius.
FOR REFERENCE: Start by aligning top of the back panel with the middle of the ear.

5. TIGHTENING

A. While holding the front panel in place, attach the loop strap on each side.
B. To tighten, anchor your fingers in on the collar and peel back the loop strap on one side, tighten and reattach. Tighten the other side equally and repeat as needed to ensure a snug, symmetrical fit.
C. Check to see that the back of the chin piece is not pressing inward on the throat. If it is, lift the back ends of the side panels up and off the trapezius muscles.

6. TIPS

A. The patient's chin should be close to the front of the chin piece but the back of the chin piece should not press on the patient's throat.
B. All slack between the front and back panels must be removed to ensure a secure fit.
C. The patient often changes position, so it may be desirable to readjust the height and the back panel straps.

INTENDED USE

To provide motion restriction of the cervical spine.

INDICATIONS

Decompressive procedures, Fracture management, Posturgical stabilization, Whiplash, Spains/strains, Cervicogenic headache, Cervical disc syndrome, Musculoskeletal pain, Radiculopathy.

CONTRAINDICATIONS

Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS

READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS

Hand wash only at 30° C. Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS

ABS, Nylon, HDPE, Polycarbonate, Deirin, Polyurethane, Cotton.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL

Product may be safely disposed in accordance with local laws.

REPORTING

Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY

Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:


- (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
- (b) The product is altered;
- (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
- (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIMITATIONS ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.


This warranty gives you specific legal rights and you may have other rights which vary from location to location.







A listing of authorized distributors in your area is available upon request.

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- (d) El producto no se utilizó de acuerdo a las Instrucciones de Uso de Aspen Medical Products, LLC.