







ASPEN HORIZON™ 456 TLSO

THORACIC LUMBAR SACRAL ORTHOSIS
LUMBAL-SAKRAL-THORAKALE ORTHESE
ORTESIS TORÁCICA LUMBOSACRALE
ORTOSI TORACO-LOMBO-SACRALE
ORTHÈSE THORACIQUE LOMBO-SACRÉE
LUMBALE SACRALE THORACALE ORTHESE
BEDERNÍ KRÍŽOVÁ HRUDNÍ ORTÉZA
LÆNDEHVIRESØJLEN THORAXORTOSE
ΟΣΦΥΪΚΗ ΘΩΡΑΚΙΚΗ ΟΡΘΩΣΗ ΙΕΡΟΥ ΟΣΤΕΟΥ
LUMBAL BRYSTORTROSE
ORTEZA PIERSIOWO-LĘDZIOWO-KRZYŻOWA
LANNERANGAN SAKRAALIN RINTAKEHÄORTOSI
LUMBAL SAKRAL THORAKAL ORTOS
ОРТЕЗ ДЛЯ ГРУДНОГО И ПОСВИЧНО-КРЕСТЦОВОГО ОТДЕЛОВ ПОЗВОНОЧНИКА
ÓRTESE TORÁCICA SACRO LOMBAR

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ENGLISH

1. SIZE ADJUSTMENT
A. To change the size of the Horizon, peel up the inside end of the side panel on each side. Adjust to proper fit and reattach.
B. The Horizon Brace sizing follows closely to standard pants sizing. Size ranges are only guidelines. Your experience and personal assessments could make the difference in choosing the correct size.

| SIZING GUIDE | | |
|--------------|---------------|--------------|
| Size | Circumference | |
| 1 | 24 – 34 in | 61 – 86 cm |
| 2 | 32 – 40 in | 81 – 102 cm |
| 3 | 38 – 46 in | 97 – 117 cm |
| 4 | 44 – 50 in | 112 – 127 cm |

2. PREPARATION
Place the pull-tabs on the belt 2 inches (5 cm) from the plastic in the back of the Horizon 456. Pull the ends of the belt away from each other ensuring that the belt is fully extended.

3. APPLICATION
A. Apply the fully extended brace making sure the back panel is centered and low on the back. Apply the brace tightly around the abdomen, the right side overlapping the left.
NOTE: When applied properly, both the upper and lower edge of the brace should be snug. If there is any gap, reapply the brace and angle the ends to insure a snug fit at both the upper and lower edge. For certain patients, the ends must be steeple angled.
B. Slip each thumb into the pull tabs. Extend the tabs out and away from the brace. When desired compression is achieved, place the hook side of the pull tabs down onto the loop material to secure.

4. ADJUSTMENT
A. Adjust the vertical height with the posterior adjustment panel. The yoke should sit at level C7 on patient.
B. Pull and tighten the shoulder straps until properly fit.
C. Fasten and tighten the chest buckle to draw the shoulder straps away from the armpits which may be more desirable for some patients.
NOTE: After initial fitting, have patient sit, then stand, and re-adjust fitting as necessary.

5. MODIFICATION OPTIONS
This product may be substantially modified by an individual with expertise. For more information, please visit www.aspenmp.com/education/product-modifications.

INTENDED USE
To provide motion restriction and stability of the Thoracic-lumbar-sacral spine.

INDICATIONS
Postoperative support, Degenerative disc disease, Bulging or herniated disc, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Kyphosis, Spondylolysis (Osteoarthritis/Degenerative Joint Disease).

CONTRAINDICATIONS
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

PRECAUTIONS
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

CARE INSTRUCTIONS
Hand wash only at 30° C. Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

MATERIALS
Nylon, Polyester, HDPE, PVC, Polyurethane, UHMWPE.

NOT MADE WITH NATURAL RUBBER LATEX.

DISPOSAL
Product may be safely disposed in accordance with local laws.

REPORTING
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

LIMITED WARRANTY
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or part(s).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:
(a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
(b) The product is altered;
(c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
(d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.
THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE, OR DEFECT IN THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.
This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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