






## ASPEN ELITE™ COMFORT

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### ENGLISH

**1. PREPARATION**  
Select appropriate brace size by measuring around the waist and hips.

SIZING GUIDE		
Size	Circumference	
X-SMALL	25.5 – 29.5 in	65 – 75 cm
SMALL	29.5 – 33.5 in	75 – 85 cm
MEDIUM	33.5 – 37.25 in	85 – 95 cm
LARGE	37.25 – 41.25 in	95 – 105 cm
X-LARGE	41.25 – 45.25 in	105 – 115 cm
XX-LARGE	45.25 – 49.25 in	115 – 125 cm

**2. ADJUST LUMBAR PAD**  
Position the lumbar pad in the desired location on the inside of the belt as shown. Press down on the hook and loop to secure.

**3. PLACEMENT**  
**A.** Apply the fully stretched brace making sure the back panel is centered and low on the back.  
**B.** Overlap the ends of the brace so the hook material grips the opposite side loop as tightly as possible.  
**C.** Be sure the top and bottom of the brace is snug. If needed, angle ends to create a snug fit.

**INTENDED USE**  
To provide motion restriction and stability of the lumbar-sacral spine.

**INDICATIONS**  
Postoperative support, Deconditioned trunk musculature, Post laminectomy syndrome, Surgical Decompression procedures, Non-surgical Decompression procedures, Degenerative disc disease, Bulging or herniated disc, Non-specific low back pain, Chronic low back pain, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Ligament sprains / Muscle strains, Spondylolysis (Osteoarthritis/Degenerative Joint Disease), Facet syndrome, SI dysfunction, Sacroiliitis, Radiculopathy Injections, Radiofrequency Ablation, Lumbar Multifidus Muscle Atrophy.

**KONTRAINDICATIONS**  
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

**PRECAUTIONS**  
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

**CARE INSTRUCTIONS**  
Hand wash only at 30° C; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

**MATERIALS**  
Nylon, Polyester, Polyethylene, Polyurethane, Acetal, PVC

**NOT MADE WITH NATURAL RUBBER LATEX.**

**DISPOSAL**  
Product may be safely disposed in accordance with local laws.

**REPORTING**  
Please report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

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## ENGLISH

### LIMITED WARRANTY

Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or part(s).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:

- (a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;
- (b) The product is altered;
- (c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product; or
- (d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OF, OR DEFECT IN, THE PRODUCT. ASPEN MEDICAL PRODUCTS, LLC SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE, AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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©2021 Aspen Medical Products, LLC. All Rights Reserved.
- (d) das Produkt nicht entsprechend der Gebrauchsanweisungen von Aspen Medical Products, LLC verwendet wird.
### ASPEN MEDICAL PRODUCTS GmbH
- (a) Produkt nie został zakupiony w firmie Aspen Medical Products, LLC ani poprzez autoryzowane kanały dystrybucji;
- (c) Do produktu wmontowano elementy, które nie zostały dostarczone przez Aspen Medical Products, LLC lub
- (d) Produkt nie jest stosowany zgodnie z instrukcjami użytkownika firmy Aspen Medical Products, LLC.